

QuickONE GEN2



QuickONE GEN2
smallest desktop
USB 2D barcode scanner

Jetzt NEU
im ALU-Gehäuse

Data Collection

- mobile 2D barcode data collector
- optional sample to patient assignment

Standard Export

- export files in csv, pdf and html format

User Interface

- Intuitive operation with touch optimized interface




QuickONE GEN2 is the smallest desktop USB barcode reader for 2D barcodes with integrated optical unit. It is possible to read single tubes with data matrix barcodes on the bottom and collect their data for later use.

The enhanced export module allows to directly generate collection reports with the captured codes and their individual timestamps for documentation.

The main parts of the QuickONE GEN2 Edition

USB desktop 2D barcode reader

The desktop 2D Barcode Reader allows scanning of barcodes directly on the desktop PC.



data collector

Captured barcodes are listed in a capturelist with their individual timestamps.


With the enhanced data collection module, it is also possible to collect additional patient information (like patient id, name...) for each sample.


export

All captured data can be exported to CSV, TXT, PDF and HTML files for documentation.

virtual keyboard

to directly insert collected barcodes in other applications, it is possible to activate the VirtualKeyboard to emulate a keyboard input in the actual active application on barcode collection.





0087656643 - Editor

QuickONE Report

file:///D:/USER/MarcoLanzinger/Desktop/...

QuickONE Export

Id: 0087656643
Name: Mayer, Markus
Birthday: Freitag, 23. Februar 1979

200000004131	DataMatrix	2014-11-14T11:48:49.2587
200000004055	DataMatrix	2014-11-14T11:48:54.2475
200000004049	DataMatrix	2014-11-14T11:48:55.0612
200000004058	DataMatrix	2014-11-14T11:48:57.6420
200000004051	DataMatrix	2014-11-14T11:48:58.7019
200000004139	DataMatrix	2014-11-14T11:49:02.4473
200000004137	DataMatrix	2014-11-14T11:49:03.1874
200000004138	DataMatrix	2014-11-14T11:49:05.7003
200000004134	DataMatrix	2014-11-14T11:49:06.5220
200000004141	DataMatrix	2014-11-14T11:49:09.0659
200000004133	DataMatrix	2014-11-14T11:49:09.7084

QuickONE Export

generated / developed by www.labo-diagnostics.de - QuickONE

ID: 0087656643
Name: Mayer, Markus
Birthday: Freitag, 23. Februar 1979

Barcode	Type	TimeStamp
200000004131	DataMatrix	2014-11-14T11:48:49.2587187+01:00
200000004055	DataMatrix	2014-11-14T11:48:54.2475956+01:00
200000004049	DataMatrix	2014-11-14T11:48:55.0612271+01:00
200000004058	DataMatrix	2014-11-14T11:48:57.642031+01:00
200000004051	DataMatrix	2014-11-14T11:48:58.7019613+01:00
200000004139	DataMatrix	2014-11-14T11:49:02.4473344+01:00
200000004137	DataMatrix	2014-11-14T11:49:03.1874019+01:00
200000004138	DataMatrix	2014-11-14T11:49:05.7003304+01:00
200000004134	DataMatrix	2014-11-14T11:49:06.5220367+01:00
200000004141	DataMatrix	2014-11-14T11:49:09.0659102+01:00
200000004133	DataMatrix	2014-11-14T11:49:09.7084032+01:00

Technical Specifications

Data Collector

Operating System	Windows 7, Windows 8, Windows 8.1, Windows 10
2D Barcode Reader	Standard and inverted 2D DotMatrix Barcodes
Export Formats	PDF, CSV, HTML
Data Exchange	USB2.0
Dimensions (total)	29 x 41 x 42 mm
Weight (total)	57 g

Order Information

QuickONE GEN2

Item-Nr.	Description
QO-00GEN2	mobile data collector for 2D barcodes with integrated optical unit

QuickONE GEN2 Options

Item-Nr.	Description
QO-PATIENT	optional patient module to capture additional patient information and data for each sample.

Our Distributor



LAB-Biotech GmbH
Eythstr. 5
89075 Ulm
Germany

Tel +49(0)731 - 1441033
Mail info@lab-biotech.de
Web www.lab-biotech.de