


one second  
**READY to SCAN !**



**DEFROST**  
your samples



### Technical Specifications

#### Defrost Roll

<b>Material</b>	3D-Printed PLA solvent resistant foam
<b>Dimensions</b>	160 x 85 x 45 mm
<b>Weight (total)</b>	100 g

Item Nr.	Description
QS-DEFROL	complete Defrost Roll
QS-REPROL	replacement foam

**1** Apply some alcohol so  
that the foam is wetted



**2** Move the rack with little  
pressure over the roll  
to defrost the rack

